COMPONENTI DEL KIT

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REGENFLEX® STARTER

Is an injectable product based on a buffered physiological solution of hyaluronic acid sodium salt (1.6%), showing notable, visco-elastic, properties highly purified, with a molecular weight between 800 and 1200 kDalton.

Main ingredient: Hyaluronic acid sodium salt.

Other ingredients

Sodium chloride, Phosphate buffer, water for injections.

Indications

REGENFLEX® STARTER is a substitutes for synovial fluid, capable of restoring the physiological and rheological properties of compromised joints, in case of pain or reduced mobility caused by degenerative or post-traumatic diseases.

Dosage and instructions for use**REGENFLEX® STARTER**

Inject 2 ml of product into the joint once a week for three weeks, subject to the physician's approval and according to health conditions of the patient.

Whenever further cycles of treatment are necessary, a six-month interval is recommended. Depending on the doctor's opinion and the patient's health conditions, more than one joint may be treated at a time.

Usage

Remove by aspiration any articular effusion before proceeding to the injection of REGENFLEX® or REGENFLEX® STARTER. Carefully extract the cap from the syringe in order to avoid any contact with the opening. Insert a proper needle twisting it on tightly so to ensure that it is airtight and to prevent leakage of the solution.

Inject the product only inside the synovial space, at room temperature and in conditions of complete asepsis.

The volume scale on the syringes is merely indicative.

REGENFLEX® BIO-PLUS

Is a buffered physiological solution, containing 2,5% of cross-linked hyaluronic acid with significant viscous and elastic properties.

The sodium hyaluronate is a fundamental component of the synovial fluid, being responsible for the particular viscous and elastic properties of the latter.


Content

1 pre-filled mono-use syringe (75 mg of reticulated hyaluronic acid in 3 ml buffered physiological solution).

Moist heat sterilized syringe.

Composition

Main ingredient: cross-linked hyaluronic acid.

Other ingredients

Sodium chloride, phosphate buffer, water for injectable solutions.

Indications

REGENFLEX® BIO-PLUS is a substitute for the synovial fluid that can restore physiological and rheological properties of compromised joints, when there are pain or reduced mobility caused by degenerative or post-traumatic diseases.

Dosage and instructions for use

Inject 3 ml of the product into the joint every 6 months, unless otherwise decided by the physician; each treatment should be tailored to the individual patient. Depending on the physician's opinion and the patient's health, it is possible to contemporarily treat several joints. The values of the volume graduation scale printed on the syringes are indicative

Mode of use

Before injecting REGENFLEX® BIO-PLUS, remove by aspiration any eventual articular effusion. Carefully extract the cap from the syringe in order to avoid any contact with the opening. Insert an appropriate needle (from 18 to 22G) twisting it on to ensure that it is air tight and to prevent leakage of the solution.

Inject the medical device only into the synovial space, at room temperature and in complete asepsis.

Warnings and contraindications REGENFLEX® STARTER and REGENFLEX® BIO-PLUS

Do not use after the expiry date printed on the box.

Do not use if the box is damaged or open.

The skin where the injection is to be carried out must be sound.

Do not carry out intra-vascular injections.

Do not inject into the knee joint when there are venous or lymphatic stasis of the limb present.

Do not inject outside the joint space, into the synovial tissue or into the joint capsule.

Do not administer in presence of excessive intra-articular effusions.

REGENFLEX® STARTER and REGENFLEX® BIO-PLUS must not be injected in case of cutaneous inflammations, infections or pathological diseases in the area to be treated.

After intra-articular injection the patient should be advised to avoid intensive physical activity and to take up normal activities only after a few days. The content of the pre-filled syringe is sterile. The external surface of the syringe is not sterile. Do not sterilize again. The device is disposable. Store at a temperature lower than 25 °C and away from heat sources. Do not freeze. Once opened the device must be used immediately and eliminated after use. Keep out of the reach of children.

Side effects

Extra-articular injections may cause local side effects. Where the needle is inserted there can be symptoms like pain, heat, reddening and swelling. In this case it is useful to apply ice to the joint treated. In normal conditions these symptoms disappear in a short time. The doctor must ensure that the patient informs him of any side effects that may eventually occur after the treatment.

In rare cases, the product can cause allergic reactions.

REGENFLEX® STARTER AND REGENFLEX® BIO-PLUS ARE ONLY SOLD UPON A MEDICAL PRESCRIPTION.

INTRA-ARTICULAR INJECTION MUST ONLY BE CARRIED OUT BY A DOCTOR. In rare cases, the product can cause allergic reactions.

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